


WHEELCHAIR BACKS



Correctly fitted lateral support, will enhance trunk stability and pelvic positioning. Our range of wheelchair backs will help reduce fatigue and maximise function for day long comfort – *work, rest or play.*

- Customisable for maximum surface support - to accommodate fixed postures such as kyphosis or lordosis
- Posterior and lateral thoracic support options - for gentle, respectful positioning
- Balancing and extending the trunk with a seat back support - will decrease pressure against the low and middle back, which in turn, decreases seat shear pressures



Slumped

Normal

Supported

Helping Hand
Pressure Care

HH a division of The Helping Hand Company



Lateral Support Back

- Compact design minimises impact on seat depth, maximising the seated area for optimum weight distribution
- Accommodates current posture with lateral support to keep body in midline
- Early preventative postural support or interim prescription

15-20" width, 12-16" high, 1.5" deep
38-50cm width, 30-40cm high, 4cm deep

Made to measure service available
Navy or black spacer cover as standard



Functional Support Back

- Promotes good 'chest open' position to facilitate body function e.g. breathing and digestion
- Adjustable lateral width to accommodate range of clients
- Continuous curved lumbar support for day long comfort

17" wide, 15" high, 4" deep
43cm wide, 38cm high, 10cm deep

Navy or black spacer cover as standard



Lowzone Back

- Designed for maximum freedom of movement, whilst providing a more functional seating posture
- Offset lumbar reduces lower back fatigue, and accommodates various heights of clients
- Adult and paediatric options

16-18" width, 16-18" high, 1 or 1.5" deep
40- 46cm width, 40-46cm high, 2.5 or 4cm deep

Black spacer, black wipedown or navy towelling covers as standard

 **COOLOVER^{TR3}** COVER OPTIONS AVAILABLE

